

Introducing the new

VisionAire™

Oxygen Concentrator

Inspired by a Glimpse into the Future of Home Care

A visionary new oxygen product has entered the scene. One that embodies all the must-have features and defines the new standard in excellence in 5 liter concentrators. Designed by AirSep® Corporation, the premier concentrator manufacturer, VisionAire™ is a brilliant technological advancement, combining the very best attributes and most desirable features.

Highlights:

- Nearly silent operation
- Significantly smaller size – both compact and lightweight
- Virtually maintenance-free
- Optional oxygen monitor
- Exceptional patient safety features and options

Clairvoyance in an Evolving Marketplace

With a clairvoyance that only a visionary manufacturer like AirSep can bring to the marketplace, this insightful concentrator design includes:

- Contemporary styling for a new generation in at-home care
- Power-Purge™ processing for greater concentrator efficiency
- Unsurpassed quality and reliability as the industry has come to expect from AirSep — the visionary leader in PSA oxygen technology.

Pure O2, Viking House, 71 Princess Road,
Urmston, Manchester, United Kingdom M41 5ST

T:- 0870 712 02 02 (UK) E: info@healthoxygen.com

W: www.healthoxygen.com



SPECIFICATIONS*

WEIGHT:

30 lb

13.6 kg

SOUND LEVEL:

40 dBA

POWER CONSUMPTION:

290 W

DIMENSIONS:

14.1 in. W x 11.5 in. D x 20.8 in. H

35.8 cm W x 29.2 cm D x 52.8 cm H

POWER:

120 VAC, 60 Hz

230 VAC, 50 Hz

ALARMS:

Power failure

High and low pressure

Temperature

Low concentration

(with optional oxygen monitor)

*Specifications subject to change without notice.



pureO₂
live longer and stronger with pure o₂