


# INOGEN ONE<sup>G3</sup> System



Inogen One® G3  
AC Power Supply\*  
Item # BA300



Inogen One® G3  
Carry Bag\*  
Item # CA300



Inogen One® G3  
Single Lithium Ion Battery\*  
Up to 4 hours run time  
Item #BA300



Inogen One® G3  
Double Lithium Ion Battery\*\*  
Up to 8 hours run time  
Item #BA316



**Only  
2.2 kg!**

## The Inogen One® G3 System (IS-300):

- 45% Smaller
- 30% Lighter
- 27% More Oxygen Per Pound  
Compared to the Inogen One G2

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere<sup>®</sup>

The new Inogen One<sup>®</sup> G3 delivers the independence of previous Inogen One portable concentrators, but it's smaller, lighter, and makes more oxygen per pound. It's a portable concentrator light enough to be worn care-free by the active oxygen user. That means there are no tanks to refill and no worry about running out of oxygen.

The Inogen One<sup>®</sup> G3 is your one solution, your single solution, for oxygen at home or away. The Inogen One<sup>®</sup> G3 is Oxygen. Anytime. Anywhere<sup>®</sup>.

**To learn more, call +31 346 216 985.**

<b>Weight</b>	2.2 kg (includes single battery)
<b>Size</b>	Length: 22 cm Width: 8 cm Height: 18 cm
<b>Oxygen Flow</b>	Intelligent Delivery Technology Four flow settings 1, 2, 3, 4
<b>Power Sensing</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power: for mobile use in car
<b>Battery</b>	Duration (single battery): 3.5 to 4 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): 7 to 8 hours Recharge approximately 4-8 hours with AC or DC power
<b>Noise Level</b>	42 Decibels*
<b>Warranty</b>	3 year**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 use, at home or away

\*At setting 2

\*\*See Inogen One<sup>®</sup> G3 Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive  
Goleta, California 93117  
1-877-845-2064  
[www.inogenone.com](http://www.inogenone.com)  
[info@inogen.net](mailto:info@inogen.net)

Distributed by:



MedSupport Systems BV  
De Corridor 12<sup>R</sup> • 3621 ZB Breukelen • NL  
[info@meds-s.com](mailto:info@meds-s.com)  
[www.medsupport-systems.com](http://www.medsupport-systems.com)

T +31 (0)346 21 69 85  
F +31 (0)346 21 66 94