

10 Litre Oxygen Concentrator

“Higher flow capabilities in a compact size”



CE

Our new 10 litre oxygen concentrator has been developed to offer optimal oxygen therapy for those requiring varying flow rates. It utilises the latest advancements in oxygen therapy and is compatible to be used with cylinder transfill systems.

www.drivedevilbiss.co.uk

drive
DeVilbiss[®]
HEALTHCARE

The compact design sits within the same footprint as our 5 litre 525 oxygen concentrator and only weighs 19kg; making this unit lightweight and compact in comparison to other 10 litre concentrators available today. Capable of delivering 87% to 96% of O₂ purity from 2 to 10 LPM, the 1025 oxygen concentrator is an ideal investment for your patients and business alike.

Features:

- » Assembled in the USA and utilises approximately 80% of the same reliable components we use for our 5 litre 525 oxygen concentrator
- » Auxiliary oxygen port – at the back of the device which can be used to fill cylinder transfill systems
- » Small, lightweight and compact design. Only weighting 19kg and within the same housing as our 5 litre 525 oxygen concentrator
- » Durable and reliable device – The unit has been tested under a number of conditions like high temperatures and high humidity
- » The unit has a bright LED panel for patient controls, with a protected cannula fitting and recessed humidifier nook to prevent damage
- » Lockable flow meter and the unit accommodates up to 15m of oxygen tubing plus 7ft of cannula



Caddy for transfill system



Auxiliary Oxygen

Patient Safety Features

Each unit comes with the DeVilbiss OSD (oxygen sensing device) which monitors the oxygen and ensures patient safety

Oxygen outlet incorporating a fire protection adapter

Visual and audible alarms built into the device indicate – power failure, high gas temperature, high pressure, low/high flow, low oxygen and service required



LED panel with hour metre

Specifications:

| | |
|-------------------------|------------------------------------|
| Outlet Pressure | 20.0 ± 1.0 psi (138 kPa +/- 7 kPa) |
| Oxygen Percentage | 2-10 LPM = 87%-96% |
| Electrical Requirements | 230V, 50 Hz |
| Power Consumption | 670 watts average |
| Weight | 19 kg |
| Dimensions | 62.2 x 34.2 x 30.4 cm |
| Warranty | 3 years |
| Classification | II Type B |

Product Code:

1025UK (10 Litre Oxygen Concentrator)

Leaflet Code: LL???

If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.